




Contents

Product	Cat. No.
Neon™ NxT Electroporation System	NEON1



Product description

The Neon™ NxT Electroporation System is a benchtop electroporation device that employs an electroporation technology which uses a pipette tip as an electroporation chamber to efficiently transfect mammalian cells including primary and immortalized hematopoietic cells, stem cells, and primary cells. The system efficiently delivers nucleic acids, proteins, and siRNA into all mammalian cell types including primary and stem cells with a high cell survival rate. The transfection is performed using as few as 1×10^4 or as many as 5×10^6 cells per reaction in a volume of 10 μL or 100 μL for a variety of cell culture formats (60 mm, 6-well, 48-well, 24-well, and 96-well).




Neon™ NxT Electroporation System

- | | |
|------------------------------|------------------------------------|
| 1. Neon™ NxT device | 4. Tube chamber for Neon™ NxT Tube |
| 2. Neon™ NxT Pipette Station | 5. Neon™ NxT Tube |
| 3. Neon™ NxT Pipette | 6. Connector cable |


For Research Use Only. Not for use in diagnostic procedures.

Instrument setup

- Place the Neon™ NxT device on a level surface. Keep the area around the unit clear for proper ventilation.
- Place the Neon™ NxT Pipette Station near the Neon™ NxT device.
- Connect the high voltage ① and low voltage ② interface connectors on the cable to the high voltage ③ and low voltage ④ interface ports at the rear of the Neon™ NxT Pipette Station.
- Connect the high voltage ① and low voltage ② interface connectors on the cable to the high voltage ③ and low voltage ④ interface ports at the rear of the Neon™ NxT Pipette Station.
- Ensure the AC power switch is in the **Off** position.
- Attach the power cord ⑤ to the AC inlet ⑥ at the rear of the Neon™ NxT device, then plug the cord into an electrical outlet.



Step 3



Step 4-6



Online resources






- Visit our product pages for protocols, safety, and additional product information.
- Go online to view related transfection products.
- For support, visit thermofisher.com/support.



Limited product warranty: Life Technologies Corporation and/or its affiliate(s) warrant their products as set forth in the Life Technologies' General Terms and Conditions of Sale found on Life Technologies' website at www.thermofisher.com/us/en/home/global/terms-and-conditions.html.

Disclaimer: TO THE EXTENT ALLOWED BY LAW, LIFE TECHNOLOGIES AND/OR ITS AFFILIATE(S) WILL NOT BE LIABLE FOR SPECIAL, INCIDENTAL, INDIRECT, PUNITIVE, MULTIPLE OR CONSEQUENTIAL DAMAGES IN CONNECTION WITH OR ARISING FROM THIS DOCUMENT, INCLUDING YOUR USE OF IT.

Manufacturer: Life Technologies Holdings Pte Ltd Block 33 | Marsiling Industrial Estate Road 3 | #07-06, Singapore 739256

©2022 Thermo Fisher Scientific Inc. All rights reserved. All trademarks are the property of Thermo Fisher Scientific and its subsidiaries unless otherwise specified.

Symbol	Action
Calculate electroporation volumes	
 Cell count calculator	a. Select Cell count calculator b. Enter the experimental parameters c. Select Calculate
Create a protocol	
 Quick run	a. Select Set up run > Quick run b. Enter the electroporation parameters c. Select Save protocol
Create a plate map	
 Create plate	a. Select Set up run > Create plate b. Select plate type (6, 12, 24, or 96-well) c. Select run order (top to bottom/by order of increasing voltage) d. Select wells and assign protocols e. Select Save plate
Edit a plate map	
 Open plate	a. Select Set up run > Open plate b. Select the protocol to modify c. Select wells and assign protocols d. Select Save plate
Open a protocol or plate map	
 Library	a. Select Library > Protocol Library or Plate Library a. Select a protocol or plate map

Symbol	Function
Open Sign in screen	
 Sign in	Select Sign in to log in to the instrument Provides access to the following items: <ul style="list-style-type: none"> ▪ Create user profile ▪ Change PIN ▪ Manage user profiles (Administrator only)
Open Settings screen	
 Settings	Select Settings Provides access to the following items: <ul style="list-style-type: none"> ▪ About instrument <ul style="list-style-type: none"> ▪ EULA ▪ Check updates ▪ Instrument settings <ul style="list-style-type: none"> ▪ Instrument name ▪ Date & time ▪ Sleep mode ▪ Brightness ▪ Network configuration ▪ Cloud region ▪ Auto sign out ▪ Maintenance & services <ul style="list-style-type: none"> ▪ Software update ▪ Self verification test ▪ Export instrument log ▪ Restore factory settings ▪ Run history <ul style="list-style-type: none"> ▪ Export run report ▪ Delete run report